


Akuport D-tec

BIORESONANCE: POTENTIZER



The Akuport D-tec produces electronic potencies from native or virtual potencies. It works just in combination with the bioresonance device Akuport MR2 or the Kindling Vistron.

The Akuport D-tec produces electronic potencies that allow the user to use a wide diversity of complementary processes. With the Akuport D-tec you can determine which potency is necessary in every individual case.



Akuport-Tower complemented with Akuport D-tec

Akuport D-tec

Bioresonance: Potentizer



- *Potency ranges: D3, D4, D6, D8, D10, D12, D15, D30, D60, D100, D200
- Output for transferring the electronic potency to a carrier medium
- Potentizing control through LED display
- * Power supply from mains with medical power supply
- *Built-in module interface for easy expansion in a plug & play-system


***Kindling Special! You will only find these features in this price range at Kindling!**

Applications

- Preparation of individual potencies in combination with bioenergetic methods (EAV, BFD, etc.)
- Preparation of auto-nosodes (e.g. with sekrete, exkrete)

Special Supplemental Options

- Combinable to Kindling Vistron
- Combinable to Akuport MR2 by plug&play-System
- Extensive special accessories, e.g. screened output beaker



Akuport D-tec + Retec-Tower



Akuport D-tec + Akuport MR2



Akuport D-tec + Vistron

Scope of delivery

In combination with Akuport MR2

Power cord
(onward transmission)
Operating instructions

In combination with Vistron or foreign devices

Medical power supply
Connection cable
Cover cap (in housing design)
Operating instructions

Dimensions (WxHxD): ca. 16x6x15cm
Weight: ca. 200g

